Panthera Dental Angulated Screw Channel Screwdriver
Instructions for use – Dental Practitioner Guide

This notice is not intended to replace or supersede sound medical judgement, dental practitioner's experience or training. Panthera Dental does not provide medical advice.

Panthera Dental Angulated Screw Channel Screwdriver
Table 1 : Screwdriver product numbers

	20mm	28mm
Unigrip	34018	34014
0.050"	34019	34015
0.048"	34020	34016
KCS	34021	34017

Description: The Panthera Dental angulated screw channel screwdriver is a reusable accessory used to install a prosthetic device with an angulated screw channel in order to restore chewing function.

Materials: The Panthera Dental angulated screw channel screwdriver is made of Phynox, a Cobalt-Chrome-Nickel alloy (ASTM 1058, ISO 5832-7, AFNOR NF S 90-403).

Package content: One (1) screwdriver.

Indications for use: The Panthera Dental screwdriver is indicated for use as an accessory to screw an angulated screw channel device on an endosseous dental implant. The Panthera Dental screwdriver is intended to be used with the Panthera Dental angulated screw channel.

Directions for use: The Panthera Dental screwdriver must be used with a torque wrench.

Table 2: Maximum angle and torque

Screwdriver	Unigrip	0.050"	0.048"	KCS
Maximum Angle	17°	20°	20°	25°
Maximum torque*	30 Ncm	15 Ncm	15 Ncm	35 Ncm

*Always use the lowest torque recommended if you have to choose between this table or the manufacturer's instructions.

It is important to apply some light pressure on the top of the screwdriver when screwing/unscrewing in order to have an optimal contact between the screwdriver head and the screw head pattern. A light pressure also has to be applied on the side of the screwdriver to align it as parallel as possible with the axis of the implant. Any original screw that was used once (screwed/unscrewed) must be replaced by a new screw. It is the dental practitioner's responsibility to verify the screwdriver head before use. If significant groove damage is observed, the screwdriver must be replaced by a new one. There is no expiry date for the screwdrivers and they can be disposed of as per local/national regulatory guidelines.

Contraindications:

- Do not screw roughly with the screwdrivers. If the limits of table 2 are not observed, stripping of the screw head pattern and the screwdriver might occur, causing irreversible damage.
- The use of a worn screw might result in damage to the screw head pattern even if all of the above recommendations are being followed.
- Using the screwdriver without proper pressure application as mentioned above (top and side) might damage the head of the screw.
- The use of a damaged screwdriver might result in damage to the head of the screwdriver or the head of the screw.
- Use only for angulated screw channels.
- Do not sterilize the screwdriver in the original plastic box. The plastic box will melt in the autoclave.

Warning: These instructions for use are not sufficient for the appropriate use of the Panthera Dental angulated screw channel screwdrivers; knowledge in dentistry and dental techniques are also required for its good use. Any use failing to comply with these instructions is considered inappropriate.

Storage: Keep inside the original box at room temperature, in a clean and dry place whenever possible.

Cleaning and sterilization: The screwdriver is sold non-sterile and is reusable. It is the responsibility of the dental practitioner to clean the screwdriver before its use, in conformance with the general dentistry practices. All Panthera Dental screwdrivers must be sterilized before use with the following validated cycle: Pre-vacuum steam sterilization for wrapped instruments, 4 minutes at 132°C, followed by a drying time of 30 minutes. The validated sterilization process must be performed in a validated sterilizer and its use should be in accordance with a recognized sterility assurance standard such as ANSI/AAMI ST79:2010/A4:2013.

Technical assistance: Our technicians are available to answer any questions on the use of Panthera Dental products. You may contact our customer service. Further information concerning the use of the angulated screw channel screwdrivers as well as the compatibility list are available on our website at www.pantheradental.com.

	Manufacturer		Keep dry
	Non-sterile		Batch code

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