## FRANÇAIS

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#### Références:

- ASTM F136-13, Standard Specification for Wrought Titanium-6Aluminum-4Vanadium ELI Extra Low Interstitial Alloy for Surgical Implant Applications (UNS F54501)
- ASTM F1472-14, Standard Specification for Wrought Titanium-6Aluminum-4Vanadium Alloy for Surgical Implant Applications
- ANSI/AAMI ST79-2010/A4-2013, Comprehensive guide to steam sterilization and sterility assurance in health care facilities
- ISO 5832-3:2016 Implants for surgery – Metallic materials – Part 3: Wrought titanium 6-aluminum 4-vanadium alloy

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**Description:** The Panthera Dental Bar is an accessory to the endosseous dental implant that supports a prosthetic device in an edentulous or partially edentulous patient, in order to restore chewing function.

**Materials:** The Panthera Dental bar is made of Ti6Al4V, grade 5, ISO 5832-3-certified titanium\*.

**Package content:** One (1) dental bar screwed onto the model, with attachments, if ordered with the bar.

**Indications for use:** The Panthera Dental Bar is indicated for use as an accessory to an endosseous dental implant to support a prosthetic device in a partially or edentulous patient for the purpose of restoring chewing function. It is intended for use to support multiple tooth prostheses in the mandible or maxilla. The prostheses may be screwed onto the implant. The Panthera Dental Milled Bars are indicated for compatibility with: Straumann S/S/P/TE RN platform Ø 4.8; S/S/P/TE WN (platform Ø 6.5); BL/BLT NC (platform Ø 3.3); BL/BLT RC platforms Ø 4.1 and Ø 4.8; the Zimmer Tapered Screw-Vent System for sizes 3.5 and 4.5, and the NobelActive: NP Ø3.5; RP Ø4.3 / Ø5.0; WP Ø5.5; NobelParallel CC: NP Ø3.75; RP Ø4.3 / Ø5.0; WP Ø5.5; NobelReplace: NP Ø3.5; RP Ø4.0 / Ø4.3 / Ø5.0; WP Ø5.0; Ø 6.0; Ø 6.0; NobelSpeedy: RP Ø4.0 / Ø5.0; WP Ø5.0 / Ø6.0; Brånemark: NP Ø3.3; RP Ø3.75 / Ø4.0; WP Ø5.0.

**Number of abutments/implants:** < 2 and > 10

**Patient's unwillingness to comply with follow-up instructions/reminders:** Patient allergic to one or more components of the attachment material

**Patient with bruxism or other uncontrolled parafunctional habit**

**Refer to the implant manufacturer's user guide for additional contraindications**

#### Side-effects:

- Irritation/inflammation for a few days
- Hypersalivation
- Temporary discomfort in the mouth
- Temporary pharyngeal reflux
- Gum sensitivity/infection

**Contraindications:** It is the responsibility of the dental practitioner to select the appropriate bar for the prosthetic project. For safe and effective use of the device, it should be installed by a trained professional. To maintain the high quality of Panthera Dental's product line, an annual bar check is recommended. A small diameter implant and an angulated abutment are not recommended for use in the posterior region of the mouth.

**Instructions for use:** Screw the bar firmly onto the implants using the appropriate torque tool. Ensure proper insertion. Refer to the screw and/or attachment manufacturer's instructions for the torque specifications required for screw and/or attachment insertion. Note that torque is never applied to the screws and attachments when shipped.

### Type 3 cobalt-based dental alloy for ceramic prostheses\* \*ISO 22674:2016 compliant

#### SCRAPE:

- Use metal bur; keep a minimum final thickness of 0.3 mm.
- Sandblast with **disposable** aluminum oxide from 110 to 250 µm with 3-bar blasting pressure.
- Steam clean and degrease with ethanol.
- Do not touch with fingers; use hemostatic forceps.

#### APPLY CERAMICS:

- It is **absolutely imperative** to use a ceramic with a thermal expansion coefficient that is as close as possible to the thermal expansion coefficient of our alloy.
- **Do not use bondor.**
- **Oxidation is not required.**
- If you decide to proceed with oxidation, place the piece in the oven for five minutes without vacuum, at the temperature of the last opaque layer.
- Apply ceramics according to the manufacturer's instructions.

#### WARNINGS AND PRECAUTIONS

In rare cases, cobalt-chromium alloys cause contact dermatitis on sensitive subjects. A contact test (patch test) is recommended.

Before installing the prosthesis, verify if there are any other metal implants in the patient's mouth. A "battery" effect can occur in the presence of different metals.

Metal dusts and fumes are hazardous to health. Use a suitable exhaust system when scraping and polishing.

We recommend that you inform your patients of the possibility that dental alloys could influence results of radiological examinations (MRT).

Table 1: PERCENTAGE CHEMICAL COMPOSITION (M/M)

COMPOSITION CHIMIQUE EN POURCENTAGE (M/M)			
Co	60	W	9
Cr	28	Si	1,5
Others Ms, Fe			
PHYSICAL AND MECHANICAL PROPERTIES			
Density	8,5 g/cm3		
Solidus-liquidus temperature	1308-1366°C		
Melting point	1420°C		
Yield load strength (Rp 0.2)	360 MPa		
Percentage elongation at fracture	16%		
Modulus of elasticity	183 GPa		
Vickers hardness	273 HV10		
Colour	White		
Thermal expansion coefficient 25-500°C	14.2-10 K -1		
Thermal expansion coefficient 25-600°C	14.4-10 K -1		
Ions release after 7 days	1,75 µg/cm <sup>2</sup>		
Maximum cooking temperature	980°C		
Suggested ceramics	Interaction Antagon by Elephant Dental BV VITA VMK Master by Vita Zahnfabrik		
Cytotoxicity test according to ISO 10993-10	Passed		
Irritation test according to ISO 10993-10	Passed		
Sensitization test according to ISO 10993-10	Passed		

### Zirconia bar - Additional information

- The metal bars of zirconia bars and bridges should be bonded only after ceramic firing and final glazing, since high temperatures destroy the composite adhesive.

- Use a good quality composite glue (Bredent type or other) and follow the manufacturer's instructions.

### Instructions for Ball Abutment - Dental Practitioner Guide

**Description:** The ball abutment is a component that screws directly into the endosseous dental implant to assist in prosthetic rehabilitation.

**Materials:** The ball abutment is made of Ti6Al4V ELI titanium, grade 23, and is ASTM F136-13\*-certified.

**Package content:** One (1) ball abutment.

**Indication:** The ball abutment is used to support a removable prosthesis in the maxilla and/or mandible. Ball abutments are compatible with the following screw-retained lines of implants. (See Table 1.)

Table 2: Panthera Dental ball abutment compatibility						
Manufacturer	System	Balls abutments Panthera Dental 3.4mm				
		platforms				
Camlog	Screw-Line Promote/Plus	3,3	3,8	4,3	5,0	6,0
	Root-Line 2 Promote/Plus	3,3	3,8	4,3	5,0	6,0
	Screw-Line Promote/Plus	3,3	3,8	4,3	5,0	
BioHorizons	Tapered Internal	4,5				
		5,7				
Nobel Biocare	NobelActive	NP				
	Replace>Select	RP				
Zimmer Biomet	Tapered SV	3,5				
	Tapered SV/AdVent	5,7				

#### Contraindications:

- Use with implants other than those identified in table 2
- Implant discrepancies > 10°
- Use for a fixed attachment
- Odd number of implants per arch
- Incorrect positioning of the implant preventing proper insertion of the prosthesis
- Refer to the implant manufacturer's User Guide for additional contraindications.

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© Panthera Dental, 2022. Tous droits réservés.