



## **PRESS RELEASE**

For immediate release

### **Panthera Dental is proud to introduce THE DIAMART™ by *Dominico Cascione, CDT.***

QUEBEC CITY, CANADA – March 2, 2016 – Panthera Dental is thrilled to announce the release of its new product: THE DIAMART™ *by Dominico Cascione, CDT.* The Diamart™ is a two-part implant bridge. The bridge, which can be milled in Zirconia, Pekkton® or PMMA, is supported by a titanium implant bar, a combination that gives the Diamart™ a long list of benefits.

The Diamart™ resolves screw channel angulation challenges. The titanium bar increases overall rigidity and offers proper support for cantilever. The easy access to the bridge increases reparability as well. This product comes with the option of a PMMA temporary and is compatible with 450+ implant platforms.

Furthermore, the restoration is 100% retrievable. The Zirconia supra structure can be segmented into multiple pieces. It provides solutions when the clinician is challenged with non-ideal implant angulations and allows a definitive restoration that satisfies the requirements of fit, esthetics and biomechanics while making the restoration retrievable.

#### **About Panthera Dental**

Panthera Dental is a world leader in CAD/CAM prosthetic implant solutions and dental sleep appliances. We are known for designing, developing, manufacturing and marketing high-level dental restoration solutions, mandibular advancement devices, and related products using superior quality materials and an advanced CAD/CAM process.

For more information, please contact our customer service at +1 855 233-0388 or visit: [www.pantheradental.com](http://www.pantheradental.com)

--30--

#### **Atts.: Images of the DIAMART**

Contact:

Karine Robichaud, Marketing Manager | +1 855 233-0388 |  
[marketing@pantheradental.com](mailto:marketing@pantheradental.com)